

Non-Contact Infrared Forehead Fever Thermometer Model GP-300

Professional Instant Digital Forehead Thermometer

Ideal for taking adult and child forehead temperatures safely and quickly with non-contact infrared technology.

Specifications

- Power supply voltage: 2 AAA batteries (included)
- Celsius/Fahrenheit availability
- Range: 89.6-109.2°F (32.0-42.9°C)
- Accuracy: +/- 0.2°C/0.4°F
- Resolution: 0.1°
- Automatic shutdown when not in use
- Size of product: 95mmx60mmx180mm
- Weight of product: 146gr
- Memory: 32 measured values
- FDA Approved

