


SafetyBLUE

PRECISION DIGITAL STEM THERMOMETERS

The Perfect Alternative to Mercury Thermometers.



Lollipop Head

Range: -58°F to 626°F and -50°C to 330°C

Resolutions: 0.1°F or °C

Accuracy: +/- 0.8°F (0.4°C) from 14 to 212°F (-10/100°C)

Auto Calibration Min/Max Memory

Auto Off Feature after 80 minutes (easily disabled for longer monitoring)

Supplied with Unique Serial Number for use as a Reference Standard

Available with an ISO 17025 Calibration Certificate

Food Safe and IP68 Waterproof

Catalog #	Stem Length (4mm O.D.)
K800-75	3" (75mm)
K800-150	6" (150mm)
K800-300	12" (300mm)
K800-460	18" (460mm)
K800-600	24" (600mm)

Bottle Thermometers (60mL)

Each is supplied with a Certificate of Accuracy



Catalog #	Application	Certified Temperature	Bottle
K1702	Block Heaters	+37°C & +95°C	Sand
K1704	Blood Banks	+4°C & +6°C	Mineral Oil
K1716	Freezers	-20°C	Mineral Oil
K1718	Refrigerators	+4°C	Mineral Oil
K1720	Incubator	+37°C	Sand
K1725	Ambient	+22°C	Mineral Oil



(631) 841-5500
 info@kesslerusa.com
 www.kesslerusa.com